

Hepatitis B Virus Surface Antigen Polyclonal Antibody PE Conjugated



Catalog No: #C91634PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hepatitis B Virus Surface Antigen Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hepatitis B virus(HBVD2)
Immunogen Type	2-80/389
Conjugates	PE
Target Name	Hepatitis B Virus Surface Antigen
Other Names	Hepatitis B Virus Surface Antigen
Accession No.	Uniprot:P03142Gene ID:944569
Uniprot	P03142
GeneID	944569
Calculated MW	31/44
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only