

B7-1/CD80 Polyclonal Antibody HRP Conjugated

Catalog No: #C91636HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	B7-1/CD80 Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC
Species Reactivity	Ms
Immunogen Type	1-100/306
Conjugates	HRP
Target Name	B7-1/CD80
Other Names	B7-1/CD80
Accession No.	Uniprot:Q00609Gene ID:12519
Uniprot	Q00609
GeneID	12519
Calculated MW	32
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only