

FITM2 Polyclonal Antibody Bio Conjugated

Catalog No: #C91638Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FITM2 Polyclonal Antibody Bio Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB
Species Reactivity	Hu Ms Rat
Conjugates	Bio
Target Name	FITM2
Other Names	FITM2
Accession No.	Uniprot:Q8N6M3Gene ID:128486
Uniprot	Q8N6M3
GeneID	128486
Calculated MW	30
Concentration	1 mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200

Note: This product is for in vitro research use only