



# Mannan Binding Lectin Polyclonal Antibody HRP Conjugated

Catalog No: #C91639HRP

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mannan Binding Lectin Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC ICC
Species Reactivity	Hu Ms Rat
Immunogen Type	1-100/248
Conjugates	HRP
Target Name	Mannan Binding Lectin
Other Names	Mannan Binding Lectin
Accession No.	Uniprot:P11226Gene ID:4153
Uniprot	P11226
GeneID	4153
Calculated MW	27
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only