

## CD4 Polyclonal Antibody HRP Conjugated

Catalog No: #C91657HRP

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD4 Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC
Species Reactivity	Hu Ms
Immunogen Type	31-130/457
Conjugates	HRP
Target Name	CD4
Other Names	CD4
Accession No.	Uniprot:P06332Gene ID:12504
Uniprot	P06332
GeneID	12504
Calculated MW	48
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:500-2000 IHC=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only