

# Mast Cell Chymase Polyclonal Antibody Biotin Conjugated

Catalog No: #C91658Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mast Cell Chymase Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB
Species Reactivity	Hu Ms Rat
Immunogen Type	28-42/247
Conjugates	Bio
Target Name	Mast Cell Chymase
Other Names	Mast Cell Chymase
Accession No.	Uniprot:P23946Gene ID:1215
Uniprot	P23946
GeneID	1215
Calculated MW	25
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200

Note: This product is for in vitro research use only