

# Mast Cell Chymase Polyclonal Antibody PE Conjugated



Catalog No: #C91658PE

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Mast Cell Chymase Polyclonal Antibody PE Conjugated                             |
| Host Species       | Rabbit  |
| Clonality          | Polyclonal  |
| Purification       | affinity purified by Protein A  |
| Applications       | Flow-Cyt  |
| Species Reactivity | Hu Ms Rat   |
| Immunogen Type     | 28-42/247   |
| Conjugates         | PE  |
| Target Name        | Mast Cell Chymase   |
| Other Names        | Mast Cell Chymase   |
| Accession No.      | Uniprot:P23946Gene ID:1215  |
| Uniprot            | P23946  |
| GeneID             | 1215  |
| Calculated MW      | 25  |
| Formulation        | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.                |
| Storage            | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

## Application Details

Flow-Cyt=1:50-200

Note: This product is for in vitro research use only