

FABP5 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91673Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FABP5 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	21-120/135
Conjugates	Bio
Target Name	FABP5
Other Names	FABP5
Accession No.	Uniprot:Q01469Gene ID:2171
Uniprot	Q01469
GeneID	2171
Calculated MW	15
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only