

## NALP3/CIAS1 Polyclonal Antibody HRP Conjugated

Catalog No: #C91674HRP

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | NALP3/CIAS1 Polyclonal Antibody HRP Conjugated                                  |
| Host Species       | Rabbit  |
| Clonality          | Polyclonal  |
| Purification       | affinity purified by Protein A  |
| Applications       | WB IHC ICC  |
| Species Reactivity | Hu Ms Rat   |
| Immunogen Type     | 501-600/1036  |
| Conjugates         | HRP   |
| Target Name        | NALP3/CIAS1   |
| Other Names        | NALP3/CIAS1   |
| Accession No.      | Uniprot:Q96P20Gene ID:114548  |
| Uniprot            | Q96P20  |
| GeneID             | 114548  |
| Calculated MW      | 114   |
| Formulation        | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.                |
| Storage            | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

## Application Details

WB=1:500-2000 IHC=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only