

Fibulin 5 Polyclonal Antibody PE Conjugated

Catalog No: #C91676PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fibulin 5 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	345-448/448
Conjugates	PE
Target Name	Fibulin 5
Other Names	Fibulin 5
Accession No.	Uniprot:Q9UBX5Gene ID:10516
Uniprot	Q9UBX5
GeneID	10516
Calculated MW	48
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only