

CPT1A Polyclonal Antibody PE Conjugated

Catalog No: #C91681PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CPT1A Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	Peptide
Immunogen Description	KLH conjugated synthetic peptide derived from human CPT1A: 1-100/773
Conjugates	PE
Target Name	CPT1A
Other Names	CPT1A
Accession No.	Uniprot:P50416Gene ID:1374
Uniprot	P50416
GeneID	1374
Calculated MW	86
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF 1:100-1:500

Note: This product is for in vitro research use only