

## NG2 Polyclonal Antibody FITC Conjugated

Catalog No: #C91683FITC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                                                 |
|--------------------|---------------------------------------------------------------------------------|
| Product Name       | NG2 Polyclonal Antibody FITC Conjugated                                         |
| Host Species       | Rabbit                                                                          |
| Clonality          | Polyclonal                                                                      |
| Purification       | affinity purified by Protein A                                                  |
| Applications       | Flow-Cyt ICC IF                                                                 |
| Species Reactivity | Hu Ms Rat                                                                       |
| Immunogen Type     | 701-800/2322                                                                    |
| Conjugates         | FITC                                                                            |
| Target Name        | NG2                                                                             |
| Other Names        | NG2                                                                             |
| Accession No.      | Uniprot:Q6UVK1Gene ID:1464                                                      |
| Uniprot            | Q6UVK1                                                                          |
| GeneID             | 1464                                                                            |
| Calculated MW      | 247                                                                             |
| Formulation        | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.                |
| Storage            | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

## Application Details

Flow-Cyt=2ug/Test ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only