

phospho-Gli1 (T101/ S102/S104) Polyclonal Antibody HRP Conjugated



Catalog No: #C91685HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	phospho-Gli1 (T101/ S102/S104) Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC ICC
Species Reactivity	Hu Ms Rat
Immunogen Type	IR(p-T)(p-S)P(p-S)SL
Conjugates	HRP
Target Name	phospho-Gli1 (T101/ S102/S104)
Other Names	phospho-Gli1 (T101/ S102/S104)
Accession No.	Uniprot:P08151Gene ID:2735
Uniprot	P08151
GeneID	2735
Calculated MW	118
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only