

## LAMP2 Polyclonal Antibody APC Conjugated

Catalog No: #C91686APC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	LAMP2 Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	301-411/411
Conjugates	APC
Target Name	LAMP2
Other Names	LAMP2
Accession No.	Uniprot:P13473Gene ID:3920
Uniprot	P13473
GeneID	3920
Calculated MW	45
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only