

Histone H4 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91687Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Histone H4 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	51-103/103
Conjugates	Bio
Target Name	Histone H4
Other Names	Histone H4
Accession No.	Uniprot:P02304Gene ID:121504
Uniprot	P02304
GeneID	121504
Calculated MW	11
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only