

CD47 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91695Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD47 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	241-323/323
Conjugates	Bio
Target Name	CD47
Other Names	CD47
Accession No.	Uniprot:Q08722Gene ID:961
Uniprot	Q08722
GeneID	961
Calculated MW	33
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only