

## TRPV1 Polyclonal Antibody Cy5 Conjugated

Catalog No: #C91698Cy5

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	TRPV1 Polyclonal Antibody Cy5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms
Immunogen Type	141-240/839
Conjugates	Cy5
Target Name	TRPV1
Other Names	TRPV1
Accession No.	Uniprot:Q8NER1Gene ID:7442
Uniprot	Q8NER1
GeneID	7442
Calculated MW	92
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only