

NEU2 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91701Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NEU2 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC IF
Species Reactivity	Ms Rat
Immunogen Type	101-200/380
Conjugates	Bio
Target Name	NEU2
Other Names	NEU2
Accession No.	Uniprot:Q9JMH3Gene ID:23956
Uniprot	Q9JMH3
GeneID	23956
Calculated MW	42
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only