

## SGOL2 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91706Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	SGOL2 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC ICC IF
Species Reactivity	Ms
Immunogen Type	401-500/1265
Conjugates	Bio
Target Name	SGOL2
Other Names	SGOL2
Accession No.	Uniprot:Q7TSY8Gene ID:68549
Uniprot	Q7TSY8
GeneID	68549
Calculated MW	145
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only