

NKG2C Polyclonal Antibody FITC Conjugated

Catalog No: #C91708FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NKG2C Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu
Immunogen Type	51-150/231
Conjugates	FITC
Target Name	NKG2C
Other Names	NKG2C
Accession No.	Uniprot:P26717Gene ID:3822
Uniprot	P26717
GeneID	3822
Calculated MW	26
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only