

DLK1 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91710Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DLK1 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB
Species Reactivity	Hu Ms Rat
Immunogen Type	251-350/383
Conjugates	Bio
Target Name	DLK1
Other Names	DLK1
Accession No.	Uniprot:P80370Gene ID:8788
Uniprot	P80370
GeneID	8788
Calculated MW	42
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200

Note: This product is for in vitro research use only