

## NAV1.7 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91712Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	NAV1.7 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	451-550/1988
Conjugates	Bio
Target Name	NAV1.7
Other Names	NAV1.7
Accession No.	Uniprot:Q15858Gene ID:6335
Uniprot	Q15858
GeneID	6335
Calculated MW	219
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200 IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only