

ENPP2 Polyclonal Antibody FITC Conjugated

Catalog No: #C91719FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ENPP2 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	601-700/915
Conjugates	FITC
Target Name	ENPP2
Other Names	ENPP2
Accession No.	Uniprot:Q13822Gene ID:5168
Uniprot	Q13822
GeneID	5168
Calculated MW	95
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only