

MyoD1/Myf3 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91722Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MyoD1/Myf3 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB
Species Reactivity	Hu Ms Rat
Immunogen Type	51-150/318
Conjugates	Bio
Target Name	MyoD1/Myf3
Other Names	MyoD1/Myf3
Accession No.	Uniprot:P10085Gene ID:17927
Uniprot	P10085
GeneID	17927
Calculated MW	35
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200

Note: This product is for in vitro research use only