

IL-4R Polyclonal Antibody FITC Conjugated

Catalog No: #C91726FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-4R Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC
Species Reactivity	Hu Ms Rat
Immunogen Type	201-300/825
Conjugates	FITC
Target Name	IL-4R
Other Names	IL-4R
Accession No.	Uniprot:P24394Gene ID:3566
Uniprot	P24394
GeneID	3566
Calculated MW	88
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only