

TNFSF14 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91728Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TNFSF14 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	141-240/240
Conjugates	Bio
Target Name	TNFSF14
Other Names	TNFSF14
Accession No.	Uniprot:O43557Gene ID:8740
Uniprot	O43557
GeneID	8740
Calculated MW	26
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only