

CD69 Polyclonal Antibody PE-Cy5.5 Conjugated

Catalog No: #C91739PE-Cy5.5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD69 Polyclonal Antibody PE-Cy5.5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	101-199/199
Conjugates	PE-Cy5.5
Target Name	CD69
Other Names	CD69
Accession No.	Uniprot:Q07108Gene ID:969
Uniprot	Q07108
GeneID	969
Calculated MW	22
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only