

ROS1 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91740Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ROS1 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	2231-2347/2347
Conjugates	Bio
Target Name	ROS1
Other Names	ROS1
Accession No.	Uniprot:P08922Gene ID:6098
Uniprot	P08922
GeneID	6098
Calculated MW	258
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only