

CD11c Polyclonal Antibody PE-Cy7 Conjugated

Catalog No: #C91743PE-Cy7

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD11c Polyclonal Antibody PE-Cy7 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	501-600/1163
Conjugates	PE-Cy7
Target Name	CD11c
Other Names	CD11c
Accession No.	Uniprot:P20702Gene ID:3687
Uniprot	P20702
GeneID	3687
Calculated MW	126
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only