

CLEC1 Polyclonal Antibody PE Conjugated

Catalog No: #C91754PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CLEC1 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	151-280/280
Conjugates	PE
Target Name	CLEC1
Other Names	CLEC1
Accession No.	Uniprot:Q8NC01Gene ID:51267
Uniprot	Q8NC01
GeneID	51267
Calculated MW	32
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only