

SLAMF7 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91755Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SLAMF7 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Ms Rat
Immunogen Type	51-160/335
Conjugates	Bio
Target Name	SLAMF7
Other Names	SLAMF7
Accession No.	Uniprot:Q9NQ25Gene ID:57823
Uniprot	Q9NQ25
GeneID	57823
Calculated MW	33
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only