

Histone H2A.x Polyclonal Antibody HRP Conjugated

Catalog No: #C91763HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Histone H2A.x Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC ICC
Species Reactivity	Hu Ms Rat
Immunogen Type	30-143/143
Conjugates	HRP
Target Name	Histone H2A.x
Other Names	Histone H2A.x
Accession No.	Uniprot:P16104Gene ID:3014
Uniprot	P16104
GeneID	3014
Calculated MW	16
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only