

## Glycophorin A Polyclonal Antibody APC Conjugated

Catalog No: #C91766APC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Glycophorin A Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Ms
Immunogen Type	74-150/150
Conjugates	APC
Target Name	Glycophorin A
Other Names	Glycophorin A
Accession No.	Uniprot:P14220Gene ID:14934
Uniprot	P14220
GeneID	14934
Calculated MW	14
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

Note: This product is for in vitro research use only