

Caspase 3 precursor Polyclonal Antibody Biotin Conjugated

Catalog No: #C91772Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Caspase 3 precursor Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	11-120/277
Conjugates	Bio
Target Name	Caspase 3 precursor
Other Names	Caspase 3 precursor
Accession No.	Uniprot:P55213Gene ID:836
Uniprot	P55213
GeneID	836
Calculated MW	28
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only