

IL15RA Polyclonal Antibody Biotin Conjugated

Catalog No: #C91776Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL15RA Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	51-150/267
Conjugates	Bio
Target Name	IL15RA
Other Names	IL15RA
Accession No.	Uniprot:Q13261Gene ID:3601
Uniprot	Q13261
GeneID	3601
Calculated MW	25
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only