

phospho-VEGFR2 (Tyr951) Polyclonal Antibody Biotin Conjugated



Catalog No: #C91797Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	phospho-VEGFR2 (Tyr951) Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	KD(p-Y)VG
Conjugates	Bio
Target Name	phospho-VEGFR2 (Tyr951)
Other Names	phospho-VEGFR2 (Tyr951)
Accession No.	Uniprot:P35968Gene ID:3791
Uniprot	P35968
GeneID	3791
Calculated MW	147
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only