

CLEC2D Polyclonal Antibody FITC Conjugated

Catalog No: #C91801FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CLEC2D Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	101-207/207
Conjugates	FITC
Target Name	CLEC2D
Other Names	CLEC2D
Accession No.	Uniprot:Q9UHP7Gene ID:29121
Uniprot	Q9UHP7
GeneID	29121
Calculated MW	10
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=2ug/Test IF=1:50-200

Note: This product is for in vitro research use only