

S100-A8 / MRP8 Polyclonal Antibody PE Conjugated

Catalog No: #C91809PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	S100-A8 / MRP8 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Ms Rat
Immunogen Type	25-89/89
Conjugates	PE
Target Name	S100-A8 / MRP8
Other Names	S100-A8 / MRP8
Accession No.	Uniprot:P05109Gene ID:6279
Uniprot	P05109
GeneID	6279
Calculated MW	10
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only