

THBS1 Polyclonal Antibody RBITC Conjugated

Catalog No: #C91815RBITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	THBS1 Polyclonal Antibody RBITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	801-900/1170
Conjugates	RBITC
Target Name	THBS1
Other Names	THBS1
Accession No.	Uniprot:P07996Gene ID:7057
Uniprot	P07996
GeneID	7057
Calculated MW	129
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only