

TLR9 Polyclonal Antibody FITC Conjugated

Catalog No: #C91816FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TLR9 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	951-1032/1032
Conjugates	FITC
Target Name	TLR9
Other Names	TLR9
Accession No.	Uniprot:Q9NR96Gene ID:54106
Uniprot	Q9NR96
GeneID	54106
Calculated MW	113
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only