

p53BP1 Polyclonal Antibody APC Conjugated

Catalog No: #C91821APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	p53BP1 Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	?ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	501-600/1977
Conjugates	APC
Target Name	p53BP1
Other Names	p53BP1
Accession No.	Uniprot:Q12888Gene ID:7158
Uniprot	Q12888
GeneID	7158
Calculated MW	213
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only