

# Acetyl Coenzyme A Carboxylase Polyclonal Antibody PE Conjugated



Catalog No: #C91823PE

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Acetyl Coenzyme A Carboxylase Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	2307-2355/2458
Conjugates	PE
Target Name	Acetyl Coenzyme A Carboxylase
Other Names	Acetyl Coenzyme A Carboxylase
Accession No.	Uniprot:O00763Gene ID:31
Uniprot	O00763
GeneID	31
Calculated MW	265
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only