

CD115/MCSF Receptor Polyclonal Antibody APC Conjugated



Catalog No: #C91826APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD115/MCSF Receptor Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	201-300/972
Conjugates	APC
Target Name	CD115/MCSF Receptor
Other Names	CD115/MCSF Receptor
Accession No.	Uniprot:P07333Gene ID:1436
Uniprot	P07333
GeneID	1436
Calculated MW	107
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only