

## EIF2AK4/GCN2 Polyclonal Antibody PE Conjugated

Catalog No: #C91829PE

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	EIF2AK4/GCN2 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	?ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1401-1500/1649
Conjugates	PE
Target Name	EIF2AK4/GCN2
Other Names	EIF2AK4/GCN2
Accession No.	Uniprot:Q9P2K8Gene ID:440275
Uniprot	Q9P2K8
GeneID	440275
Calculated MW	189
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only