

Urocortin 3 Polyclonal Antibody FITC Conjugated

Catalog No: #C91834FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Urocortin 3 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	120-157/161
Conjugates	FITC
Target Name	Urocortin 3
Other Names	Urocortin 3
Accession No.	Uniprot:Q969E3Gene ID:114131
Uniprot	Q969E3
GeneID	114131
Calculated MW	44304
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only