

RGS5 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91837Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RGS5 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	81-181/181
Conjugates	Bio
Target Name	RGS5
Other Names	RGS5
Accession No.	Uniprot:O08850Gene ID:19737
Uniprot	O08850
GeneID	19737
Calculated MW	22
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only