

P38 MAPK/ MAPK14 Polyclonal Antibody PE Conjugated

Catalog No: #C91840PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	P38 MAPK/ MAPK14 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	?ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	121-220/360
Conjugates	PE
Target Name	P38 MAPK/ MAPK14
Other Names	P38 MAPK/ MAPK14
Accession No.	Uniprot:Q16539Gene ID:1432
Uniprot	Q16539
GeneID	1432
Calculated MW	41
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only