

## CD40/TNFRSF5 Polyclonal Antibody APC Conjugated

Catalog No: #C91854APC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD40/TNFRSF5 Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	201-277/277
Conjugates	APC
Target Name	CD40/TNFRSF5
Other Names	CD40/TNFRSF5
Accession No.	Uniprot:P25942Gene ID:958
Uniprot	P25942
GeneID	958
Calculated MW	30
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

Note: This product is for in vitro research use only