

CD40/TNFRSF5 Polyclonal Antibody PE-Cy7 Conjugated



Catalog No: #C91854PE-Cy7

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD40/TNFRSF5 Polyclonal Antibody PE-Cy7 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	201-277/277
Conjugates	PE-Cy7
Target Name	CD40/TNFRSF5
Other Names	CD40/TNFRSF5
Accession No.	Uniprot:P25942Gene ID:958
Uniprot	P25942
GeneID	958
Calculated MW	30
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only