

Clathrin heavy chain Polyclonal Antibody RBITC Conjugated

Catalog No: #C91856RBITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Clathrin heavy chain Polyclonal Antibody RBITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	401-500/1675
Conjugates	RBITC
Target Name	Clathrin heavy chain
Other Names	Clathrin heavy chain
Accession No.	Uniprot:Q00610Gene ID:1213
Uniprot	Q00610
GeneID	1213
Calculated MW	192
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only